SRA HUMAN SCRUBBER





0

This resource was made possible through funding provided under the Epidemiology and Laboratory Capacity for Prevention and Control of Emerging Infectious Diseases (ELC) Cooperative Agreement (CK24-0002), Project D: Advanced Molecular Detection to the Florida Department of Health. The conclusions, findings, and opinions expressed by authors do not necessarily reflect the official position of the U.S. Department of Health and Human Services, the Public Health Service, or the Centers for Disease Control and Prevention.

OVERVIEW

Purpose

Identify and remove or mask human-origin reads from FASTQ files

Usage

 Used by public health labs and researchers to clean sequencing data by removing human-derived reads, helping ensure privacy compliance before sharing or analyzing pathogen genomes to study outbreaks and viral evolution.


Dependencies

Unix Tools, Curl and Pre-Built k-mer reference database (init_db.sh)



REVIEW OF K-MERS

- Short, fixed-length subsequences of DNA sequences
 - "K" refers to the length of the sequence (Ex. 3-mer, 4-mer, 32-mer
- Human Scrubber uses a 32-mer
 - Human scrubber uses a reference k-mer database to map query reads to known pathogen genomes





WORKFLOW





+

APPLICATION

Objective:

Use a single end fastq file and make a clean dataset using Human Scrubber

APPLICATION CONT.


cd /blue/bphl-<state>/<user>/repos/bphlmolecular/

git clone https://github.com/ncbi/srahuman-scrubber

mkdir analysis/

cd analysis/

cp /blue/bphl-<state>/<user>/repos/bphlmolecular/human-scrubber/*






+


APPLICATION CONT.

```
[n.yengalareddy@login8 sra-human-scrubber0617]$ ./scripts/scrub.sh test
2025-06-19 22:37:11
                       aligns to version 0.801
                       hardware threads: 64, omp threads: 64
2025-06-19 22:37:11
                       loading time (sec) 6
2025-06-19 22:37:17
                       /tmp/tmp.bjRmMbNeAI/temp.fasta
2025-06-19 22:37:17
                       FastaReader
2025-06-19 22:37:17
2025-06-19 22:37:17
                       100% processed
2025-06-19 22:37:17
                       total spot count: 2
2025-06-19 22:37:17
                       total read count: 2
2025-06-19 22:37:17
                       total time (sec) 6
  spot(s) masked or removed.
DB version is 20250325v2
test succeeded
```



APPLICATION CONT.


./scripts/scrub.sh -p 8 /blue/bphlflorida/n.yengalareddy/repos/bphl-molecular/analysis/srahuman-scrubber0617/test/scrubber_test.fastq






C

APPLICATION CONT.



@test.2 2 length=286
GGCAATGCATTTACAGTACAAAAGACATACTGTTCTAATGTTGAATTCACTTTGAATTTATCAAAACACTCTACACGAGCACGTGCAGGTATAATTCT
+test.2 2 length=286

A?AA1C1@1BFFDGB1GGGGGGFHCGBGHHH3DB3F33FHFGHGFHGHBFHHHHGHHHHHHHHEGBEFGHHHHHHHE?/EEAE/B/CEFGFGHDFFGHH @test.216524 216524 length=194



PULL HUMAN SCRUBBER-DOCKER

- DockerHub
 - https://hub.docker.com/r/ncbi/sra-human-scrubber
 - docker pull ncbi/sra-human-scrubber

```
Here the command is given the path to your local fastq file as argument docker run -it -v $PWD:$PWD:rw -w $PWD ncbi/sra-human-scrubber:latest /opt/scrubber/scripts/scrub.sh path-to-fastq-file/filename.fastq

Example: docker run -it -v $PWD:$PWD:rw -w $PWD ncbi/sra-human-
```


Example: docker run -it -v \$PWD:\$PWD:rw -w \$PWD ncbi/sra-human-scrubber:latest /opt/scrubber/scripts/scrub.sh MyFastqFile.fastq

```
2022-09-06 21:35:04
                       aligns_to version 0.707
2022-09-06 21:35:04
                      hardware threads: 8, omp threads: 8
2022-09-06 21:35:04
                      loading time (sec) 0
                       /tmp/tmp.Ccqruccyoq/temp.fasta
2022-09-06 21:35:04
2022-09-06 21:35:04
                       FastaReader
2022-09-06 21:35:04
                       0% processed
2022-09-06 21:35:06
                       100% processed
2022-09-06 21:35:06
                       total spot count: 216859
2022-09-06 21:35:06
                       total read count: 216859
2022-09-06 21:35:06
                      total time (sec) 2
129 spot(s) masked or removed.
```




PULL HUMAN SCRUBBER-BIOCONDA

- Bioconda
 - https://anaconda.org/bioconda/sra-human-scrubber
 - conda install bioconda::sra-human-scrubber
- On HPG:
 - module load sra_human_scrubber
 - Run *module spider sra-human-scrubber* to see what environment modules are available for sra human scrubber
- Note: SRA Human Scrubber is included in Bactopia



APPLYING HUMAN SCRUBBER

- To retroactively apply SRA Human Scrubber to your SRA submissions, email the SRA Help Desk
 - <u>sra@ncbi.nlm.nih.gov</u>
 - Request HRRT be activated for your BioProject
 - Include your BioProject Number
 - Depending on the number of samples it'll take about a week for Human Scrubber to be applies.
- Will also be applied to future submissions
 - Better to do it in-house before SRA submission to protect possible PHI breach




NCBI DATASETS



NCBI DATASETS OVERVIEW

- NCBI has data sets for almost anything you could ever want
- Taxonomy, gene, and genome level
 - Special data set for viruses
- Can access in 3 ways
 - CLI
 - GitHub
 - NCBI Website
- Excellent How-To Guides on their website
 - https://www.ncbi.nlm.nih.gov/datasets/d ocs/v2/how-tos/





NCBI DATASETS GUI

Bacteria / Pseudomonadota / Gammaproteobacteria / Pseudomonadales / Pseudomonadaceae /

Pseudomonas aeruginosa 🜣

Pseudomonas aeruginosa is a species of g-proteobacteria in the family Pseudomonadaceae.

| Browse taxonomy |
|-----------------|
|-----------------|

| NCBI Taxonomy ID | 287 | |
|-------------------------|--|--|
| Taxonomic rank | species | |
| Current scientific name | Pseudomonas aeruginosa (Schroeter 1872) Migula 1900 (Approved Lists 1980) NOMEN APPROBBATUM Type Material | |
| Basionym | "Bacterium aeruginosum" Schroeter 1872 | |

Genome

Browse all 31,553 genomes

View taxonomic details


| Nucleotide | | Protein | | |
|--------------------------|-----------|-----------------------------|------------|--|
| All nucleotide sequences | 3,148,747 | Protein sequences | 24,373,600 | |
| Genomic sequences | 3,148,394 | Conserved domains | 9 | |
| mRNA sequences | 100 | 3D structures | 2,943 | |
| GEO Datasets | | Sequence Read Archive (SRA) | | |
| Datasets | 24 | All SRA experiments | 62,293 | |
| Series | 496 | DNA | 55,759 | |
| Samples | 6,728 | RNA | 6,350 | |
| Platforms | 103 | | | |
| PopSet | | Projects and samples | | |
| Phylogenetic studies | 468 | BioProject | 3,277 | |
| Population studies | 241 | BioSample | 75,992 | |






NCBI DATASETS CLI

- Two CLI tools
 - Datasets: download sequence data across all domains of life
 - Dataformat: convert metadata from JSON to other formats
- Commands follow a standard syntax








NCBI DATASETS CLI CONT.

- 3-Step Conda install (includes both datasets and dataformats in the conda package)
 - 1. Create the conda environment: conda create –n ncbi_datasets
 - 2. Activate the environment: conda activate ncbi_datasets
 - Install the datasets conda package: conda install -c conda-forge ncbidatasets-cli
 - Note the switch from _ to in ncbi-datasets
 - Example code
 - datasets download genome accession GCA_020809405.1
 - datasets download genome accession GCA_020809405.1 GCA_020748185.1
 - An example of multiple genomes



NCBI DATASETS GITHUB

- Request a new feature or submit a bug report
 - .github/ISSUE_TEMPLATE
 - bug_report.md
 - feature_request.md
 - Current version as of 6/23/25 is v18.x



CONCLUSION





Fundamentals of Human Scrubber and NCBI Datasets



Installation and setup of Human Scrubber in HPG



Successfully executed job query for Human Scrubber



Generated output file for Human Scrubber





Advanced Molecular Detection Southeast Region Bioinformatics

Questions?

<u>bphl-sebioinformatics@flhealth.gov</u>

Molly Mitchell, PhD
Bioinformatics Supervisor
Molly.Mitchell@flhealth.gov

Nikhil Reddy, MS
Bioinformatician
Nikhil.Yengala@flhealth.gov

Sam Bernhoft, MPH
Bioinformatician
Samantha.bernhoft@flhealth.gov