



Advanced Molecular Detection

Southeast Region Bioinformatics

Outline



Updates



Agenda



Bactopia



ECTyper



EmmTyper



Questions

Updates – ABiL Trainings

Two upcoming trainings available to you

1. In-person training – tentatively scheduled for October 16-19 – thoughts?
 1. 4 days at Georgia Tech
 2. Intermediate to advanced bioinformatics – you'll need good fundamentals to attend.
2. Online Trainings – courses are still being developed; we will notify once they are available

Agenda

August 21 – Bactopia Tools: FastANI and GAMMA

September 4 **rescheduled to** **September 11** – Bactopia Tools: Hicap and HpsuisSero

Future Trainings

- ONT & FL's Flisochar pipeline
- StaPH-B Toolkit Programs/Pipelines
- GISAID flagged SARS-CoV-2
- R Training Series
- Dryad pipeline
- ...and more

Bactopia

- Bactopia is a flexible pipeline for complete analysis of bacterial genomes
- Bactopia was inspired by Staphopia, a workflow that targets *Staphylococcus aureus* genomes
- Bactopia was developed from scratch prioritizing usability, portability, and speed

Bactopia Usage

- Bactopia uses Nextflow to manage the workflow – which supports many types of environments (e.g., cluster or cloud)
- Bactopia allows for the usage of many public datasets as well as your own datasets to further enhance the analysis of your sequencing data
- Bactopia only uses software packages available from Bioconda (or other Anaconda channels) to make installation simple for *all* users



- ⚠ *Aborting poor quality samples prevents downstream failures which would stop all samples*
 - Too few reads or basepairs
 - Coverage below minimum
 - Paired-end with different read counts
 - Paired-end with skewed proportions
 - Genome size below minimum
 - Genome size exceeds maximum
 - 0 assembled contigs
 - Assembled size below minimum

Legend

- Process uses F
- Process uses C
- Process uses M
- Process uses C
- Minimum QC ne

Supplemented By Ba

- Generic datase
- Species-specifi

Bactopia Processes

- Gather Samples**
Collect local files and/or download
- QC Reads**
Trim and filter low quality reads, check coverage, and generate quality
- Minmer Sketch and Minmer**
Create minmer sketches and query against GenBank
- Call Variants**
Determine SNPs and InDels against reference
- Ariba Analysis**
Query FASTQs against Ariba database
- Mapping Query**
Align to a reference and determine
- Assemble Genome & Assemble**
Create a de novo assembly and assess the quality of the assembly
- Annotate Genome**
Predict genes and proteins from assembly
- Antimicrobial Resistance**
Identify presence of AMR and virulence
- Blast**
Align genes, proteins, or primers against
- Sequence Type**
Determine sequence type based on

Workflow

Bactopia Tools

ECTyper (an Easy Typer)

- ECTyper is a standalone versatile serotyping module for *Escherichia coli*
- Supports both *.fasta* (assembled) and *.fastq* (raw reads) file formats
- This tool provides convenient species identification coupled with quality control module giving a complete, transparent, and reference laboratory suitable report on *E. coli* serotyping

[phac-nml/ecoli_serotyping: In silico prediction of E. coli serotype \(github.com\)](https://phac-nml/ecoli_serotyping:In%20silico%20prediction%20of%20E.%20coli%20serotype%20(github.com))

Installation

Can be installed through conda

```
conda create -yp /blue/bphl-<state>/<user>/conda_envs/ectyper/  
conda activate /blue/bphl-<state>/<user>/conda_envs/ectyper/  
conda install -c conda-forge -c bioconda ectyper
```

Usage

```
thsalikilakshmi@login1:/blue/bphl-florida/thسالikilakshmi/data/HAI/20220727_jax_220708_PLN_WLK_MS_test/assemblies
(/blue/bphl-florida/thسالikilakshmi/training/conda_envs/ectyper) [thسالikilakshmi@login1 assemblies]
$ ectyper -h
usage: ectyper [-h] [-V] -i INPUT [-c CORES] [-opid PERCENTIDENTITYOTYPE]
              [-hpid PERCENTIDENTITYHTYPE] [-opcov PERCENTCOVERAGEOTYPE]
              [-hpcov PERCENTCOVERAGEHTYPE] [--verify] [-o OUTPUT]
              [-r REFSEQ] [-s] [--debug] [--dbpath DBPATH]

ectyper v1.0.0 database v1.0 Prediction of Escherichia coli serotype from raw
reads or assembled genome sequences. The default settings are recommended.

optional arguments:
  -h, --help                show this help message and exit
  -V, --version              show program's version number and exit
  -i INPUT, --input INPUT    Location of E. coli genome file(s). Can be a single
                              file, a comma-separated list of files, or a directory
  -c CORES, --cores CORES    The number of cores to run ectyper with
  -opid PERCENTIDENTITYOTYPE, --percentIdentityOtype PERCENTIDENTITYOTYPE
                              Percent identity required for an O antigen allele
                              match [default 90]
  -hpid PERCENTIDENTITYHTYPE, --percentIdentityHtype PERCENTIDENTITYHTYPE
                              Percent identity required for an H antigen allele
                              match [default 95]
  -opcov PERCENTCOVERAGEOTYPE, --percentCoverageOtype PERCENTCOVERAGEOTYPE
                              Minumum percent coverage required for an O antigen
                              allele match [default 95]
  -hpcov PERCENTCOVERAGEHTYPE, --percentCoverageHtype PERCENTCOVERAGEHTYPE
                              Minumum percent coverage required for an H antigen
```



Input File Format

Takes *.fasta* files as input

```
$ ectyper -i JBE22000155.fasta -o results_ectyper
```

Output

Ectyper serotyping results are available in a tab-delimited **output.tsv** file consisting of 16 columns:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Name	Species	O-type	H-type	Serotype	QC	Evidence	GeneScore	AlleleKey	GeneIden	GeneCove	GeneCont	GeneRang	GeneLeng	Database	Warnings
2	JBE220001	-	-	H8	-:H8	-	Based on	fliC:1;	H8-1-fliC-	100;	100;	73;	15181-166	1479;	v1.0 (11-0)	-
3									.							

Interpreting Results

- Name – Sample name (usually a unique identifier)
- Species: the species column provides valuable species identification information in case of inadvertent sample contamination or mislabeling events
- O-type: O antigen
- H-type: H antigen
- Serotype: Predicted O and H antigen(s)
- QC: The Quality Control value summarizing the overall quality of prediction
- Evidence: How many alleles in total used to both call O and H antigens
- GeneScores: ECTyper O and H antigen gene scores in 0 to 1 range
- AllelesKeys: Best matching ECTyper database allele keys used to call the serotype

Interpreting Results

- GenIdentities(%): %identity values of the query alleles
- GeneCoverages(%): %coverage values of the query alleles
- GeneContigNames: the contig names where the query alleles were found
- GeneRanges: genomic coordinates of the query alleles
- GeneLengths: allele lengths of the query alleles
- Database: database release version and date
- Warnings: any additional warnings linked to the quality control status or any other error message(s)

EmmTyper

- Emm Automatic Isolate Labeller
- emmtyper is a command line tool for emm-typing of *Streptococcus pyogenes* using a de novo or complete assembly
- Tool has two basic modes:
 - blast: contigs are blast against the trimmed FASTA database curated by the CDC
 - pcr: in-silico PCR is done on the contigs using the isPCR tool

[MDU-PHL/emmtyper: emm Automatic Isolate Labeller \(github.com\)](https://github.com/MDU-PHL/emmtyper)

How emm genes work?

- The difficulty with performing M-typing is that there is a single gene of interest (emm), and two other homologous genes (enn and mrp), often referred to emm-like
- Homologous genes may or may not occur in the isolate of interest
- When performing emm-typing from an assembly, we can distinguish between one or more clusters of matches on the contigs
- The best match for each of the clusters identified is then parsed from the BLAST results

Installation

Can be installed through conda

```
conda create -yp /blue/bphl-<state>/<user>/conda_envs/emmtyper/  
conda activate /blue/bphl-<state>/<user>/conda_envs/emmtyper/  
conda install -c conda-forge -c bioconda emmtyper
```

Usage

emmtyper has two workflows:

1. Directly BLASTing the contigs against the database
2. Using isPcr to generate an *in-silico* PCR product which is BLAST against the database

Help Menu

```
thsalikilakshmi@login1:/blue/bphl-florida/thسالikilakshmi/data/HAI/20220727_jax_220708_PLN_WLK_MS_test/assemblies
(/blue/bphl-florida/thسالikilakshmi/training/conda_envs/emmtyper) [thسالikilakshmi@login1 assemblies]
es]$ emmtyper --help
Usage: emmtyper [OPTIONS] [FASTA]...

Welcome to emmtyper.

Usage:

emmtyper *.fasta

Options:
--version                Show the version and exit.
-w, --workflow [blast|pcr] Choose workflow [default: blast]
-d, --blast_db TEXT      Path to EMM BLAST DB [default: /blue/bphl-f
lorida/thسالikilakshmi/training/conda_envs/e
mmtyper/lib/python3.11/site-
packages/emmtyper/db/emm.fna]
-k, --keep                Keep BLAST and isPcr output files.
-d, --cluster-distance INTEGER Distance between cluster of matches to
consider as different clusters. [default:
500]
-o, --output TEXT         Output stream. Path to file for output to a
file. [default: stdout]
-f, --output-format [short|verbose|visual] Output format.
--dust [yes|no|level window linker] [BLAST] Filter query sequence with DUST.
[default: no]
--percent-identity INTEGER [BLAST] Minimal percent identity of
```



Input & Results

- Takes *.fasta* files as input

```
$ emmtyper JBE*.fasta > results_emmtyper
```

- Results

	A	B	C	D	E	
1	Isolate name	No.of clusters	Predicted emm-type	emm-like genes	emm cluster	
2	GA230457.tmp	2	EMM89.0	EMM203.4	E4	

Interpreting Results

- emmtyper has three different result formats: short, verbose, and visual
- emmtyper by default produces the short version. This consists of five values in tab-separated format printed to stdout. These values are:
 - Isolate name
 - Number of clusters: should be between 1 and 3, larger values could indicate contamination
 - Predicted emm-type
 - Possible emm-like alleles (semi colon separated list)
 - EMM cluster: Functional grouping of EMM types into 48 clusters



Advanced Molecular Detection Southeast Region Bioinformatics

Questions?

bphl-sebioinformatics@flhealth.gov

Lakshmi Thsaliki, MS

Bioinformatician

Lakshmi.Thsaliki@flhealth.gov

Molly Mitchell, PhD

Bioinformatician

Molly.Mitchell@flhealth.gov