



# **Advanced Molecular Detection**

## **Southeast Region Bioinformatics**

# Outline



Agenda



Updates



Shigatyper



Shigeifinder



Questions

# Agenda

**January 08** – Bactopia Tools: spatyper and ssuissero

**January 22** – Bactopia Pipeline

## Future Trainings

- ONT & FL's Flisochar pipeline
- StaPH-B Toolkit Programs/Pipelines
- GISAID flagged SARS-CoV-2
- R Training Series
- Dryad pipeline
- ...and more

# Updates

- AMD Southeast Region Annual Needs Assessment Survey is sent out last Tuesday. Please make sure you finish the survey before **December 22<sup>nd</sup>**
- Feel free to forward the survey to any new staff member not included in this email
- Also, please have at least one staff member respond back to make any updates to your lab AMD contact list

# Shigatyper

- ShigaTyper is a quick and easy tool designed to determine *Shigella* serotypes using Illumina (single or paired-end), or Oxford Nanopore reads with low computation requirement
- ShigaTyper supports compressed FASTQs as inputs
  - FASTQs can be single-end or paired-end Illumina reads, or reads from Oxford Nanopore

[CFSAN-Biostatistics/shigatyper: CFSAN Shigella Typing Pipeline \(github.com\)](https://github.com/CFSAN-Biostatistics/shigatyper)

# Installation

Can be installed through conda

```
conda create -yp /blue/bphl-<state>/<user>/conda_envs/shigatyper/  
conda activate /blue/bphl-<state>/<user>/conda_envs/shigatyper/  
conda install -c conda-forge -c bioconda shigatyper
```

# Usage

```
usage: shigatyper [-h] [--R1 FASTA] [--R2 FASTA] [--SE FASTA] [--ont] [-n SAMPLE_NAME] [-o OUTDIR] [--verbose] [--version]
```

ShigaTyper v. 2.0.5, 2022

A WGS-based genoserotyping pipeline for *Shigella* spp.

Yun Wu, Henry K Lau, Teresa Lee, David K Lau, Justin Payne

The bacteria *Shigella* spp., consisting of 4 species and >50 serotypes, cause shigellosis, a foodborne disease of significant morbidity, mortality, and economic loss worldwide. Classical *Shigella* identification based on selective media and serology is tedious, time-consuming, expensive, and not always accurate. Molecular diagnostic assay does not distinguish *Shigella* at species level or from enteroinvasive *Escherichia coli* (EIEC). We inspected the whole genome sequencing (WGS) data from 219 *Shigella* isolates and observed low concordance rate between conventional designation and molecular serotyping, 86.8% and 78.9% at species and serotype level, respectively. Serotype determinants for 6 additional serotypes were identified. Examination of differentiation gene markers commonly perceived as characteristic hallmarks in *Shigella* showed high variability among different serotypes. Using this information, we developed ShigaTyper, an automated workflow that utilizes limited computational resources to accurately and rapidly determine 58 *Shigella* serotypes using Illumina paired end WGS reads. *Shigella* serotype determinants and species-specific



# Input

Takes paired-end reads, single-end reads, or Oxford Nanopore reads as input

```
shigatyper --SE JBI22000647.fasta > shigatyper_results
```



# Results

	Hit	Number of reads	Length Covered	reference length	% covered	Number of variants	% accuracy
0	ipaH_c	1	717	780	91.9	0	100
1	Sf_wzx	1	0	1257	0	0	
2	Sf_wzy	1	0	1149	0	0	
3	gtrl	1	0	1521	0	0	
4	Oac1b	1	0	1002	0	0	

# ShigEifinder

- This tool differentiates *Shigella*/EIEC using cluster-specific genes and identify the serotype using O-antigen/H-antigen genes
- This pipeline can serotype over 59 *Shigella* and 22 EIEC using either assembled genomes or WGS reads

[LanLab/ShigEiFinder: Cluster informed Shigella and EIEC serotyping tool from Illumina reads and assemblies \(github.com\)](#)

# Installation

Can be installed through conda

```
conda create -yp /blue/bphl-  
<state>/<user>/training/conda_envs/shigeifinder/  
conda activate /blue/bphl-<state>/<user>/training/conda_envs/shigeifinder/  
conda install -c conda-forge -c bioconda shigeifinder
```

# Usage

```
usage:
Assembly fasta input/s:
  ShigeiFinder.py -i <input_data1> <input_data2> ... OR
  ShigeiFinder.py -i <directory/*>
Paired end raw read fastq(.gz) input/s:
  ShigeiFinder.py -r -i <Read1> <Read2> OR
  ShigeiFinder.py -r -i <directory/*>
Single end raw read fastq(.gz) input/s:
  ShigeiFinder.py -r --single_end -i <Reads> OR
  ShigeiFinder.py -r --single_end -i <directory/*>

options:
  -h, --help            show this help message and exit
  -i I [I ...]          <string>: path/to/input_data
  -r                    Add flag if file is raw reads.
  -t T                  number of threads. Default 4.
  --single_end          Add flag if raw reads are single end rather than paired.
  --hits               To show the blast/alignment hits
  --dratio             To show the depth ratios of cluster-specific genes to House Keeping genes
  --update_db          Add flag if you added new sequences to genes database.
  --output OUTPUT      output file to write to (if not used writes to stdout)
  --check              To show the blast/alignment hits
  --o_depth O_DEPTH    When using reads as input the minimum depth percentage relative to genome average for positive O antigen gene call (default 1.0).
  --ipaH_depth IPA_H_DEPTH
                        When using reads as input the minimum depth percentage relative to genome average for positive ipaH gene call (default 1.0).
  --depth DEPTH        When using reads as input the minimum read depth for non ipaH/Oantigen gene to be called (default 10.0).
  --tmpdir TMPDIR      temporary folder to use for intermediate files
  --noheader           do not print output header
  -v, --version        Print version information.
```



# Input

Takes fasta files as input

```
$ shigeifinder -i JBI22000648.fasta > shigeifinder_results
```

# Results

	#SAMPLE	ipaH	VIRULENCE	CLUSTER	SEROTYPE	O_ANTIGEN	H_ANTIGEN	NOTES
1								
2	JB122000648	+	38	C3	SF1b	SF1-5		



# **Advanced Molecular Detection Southeast Region Bioinformatics**

## **Questions?**

[bphl-sebioinformatics@flhealth.gov](mailto:bphl-sebioinformatics@flhealth.gov)

**Lakshmi Thsaliki, MS**

Bioinformatician

[Lakshmi.Thsaliki@flhealth.gov](mailto:Lakshmi.Thsaliki@flhealth.gov)

**Molly Mitchell, PhD**

Bioinformatician

[Molly.Mitchell@flhealth.gov](mailto:Molly.Mitchell@flhealth.gov)