

13 Oct 2025

# Aquascope

Advanced Molecular Detection  
Southeast Region Bioinformatics

This resource was made possible through funding provided under the Epidemiology and Laboratory Capacity for Prevention and Control of Emerging Infectious Diseases (ELC) Cooperative Agreement (CK24-0002), Project D: Advanced Molecular Detection to the Florida Department of Health. The conclusions, findings, and opinions expressed by authors do not necessarily reflect the official position of the U.S. Department of Health and Human Services, the Public Health Service, or the Centers for Disease Control and Prevention.

# Overview

## Purpose:

Pipeline is developed for *SARS-CoV-2* wastewater sequencing analysis, enabling standardized qc, read trimming, alignment, primer trimming and variant deconvolution to estimate lineage abundance in wastewater samples

## Usage:

Can be used by public health labs and researchers for surveillance of *SARS-CoV-2* variants through wastewater for molecular typing, outbreak investigation, surveillance, and comparative genomic studies

## Dependencies:

Nextflow


Apptainer

Conda

Primer Scheme (ARTIC V3) – Required for Illumina Runs



# Workflow



# Application

## Objective

Use SARS-CoV-2 Illumina FASTQs paired sequences (acquired from wastewater) from NCBI database, and use Aquascope to analyze the samples

# Application Cont.

```
cd blue/bphl-<state>/<user>/repos/bphl-molecular/
```

```
git clone https://github.com/CDCgov/aquascope
```

```
mkdir analysis/
```

```
cd analysis/
```

```
cp /blue/bphl-<state>/<user>/repos/bphl-molecular/aquascope/
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/
Name
..
assets
bin
conf
docs
input
misc
modules
results
subworkflows
test
workflows
aquascope.14354509.err
aquascope.14354509.out
CHANGELOG.md
CITATIONS.md
code-of-conduct.md
CONTRIBUTING.md
DISCLAIMER.md
LICENSE
main.nf
mkdocs.yml
modules.json
nextflow.config
nextflow_schema.json
open_practices.md
README.md
ro-crate-metadata.json
rules_of_behavior.md
run_workflow.sh
tower.yml
VERSION
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/
Name
samplesheet.csv
SRR12618483_R1.fastq.gz
SRR12618483_R2.fastq.gz
SRR12618484_R1.fastq.gz
SRR12618484_R2.fastq.gz
SRR12618485_R1.fastq.gz
SRR12618485_R2.fastq.gz
SRR12618486_R1.fastq.gz
SRR12618486_R2.fastq.gz
SRR12618487_R1.fastq.gz
SRR12618487_R2.fastq.gz
SRR12618488_R1.fastq.gz
SRR12618488_R2.fastq.gz
SRR12618489_R1.fastq.gz
SRR12618489_R2.fastq.gz
SRR12618490_R1.fastq.gz
SRR12618490_R2.fastq.gz
SRR12618491_R1.fastq.gz
SRR12618491_R2.fastq.gz
SRR12618492_R1.fastq.gz
SRR12618492_R2.fastq.gz
SRR12618493_R1.fastq.gz
SRR12618493_R2.fastq.gz
SRR12618494_R1.fastq.gz
SRR12618494_R2.fastq.gz
SRR12618495_R1.fastq.gz
SRR12618495_R2.fastq.gz
SRR12618496_R1.fastq.gz
SRR12618496_R2.fastq.gz
SRR12618497_R1.fastq.gz
SRR12618497_R2.fastq.gz
SRR12618498_R1.fastq.gz
SRR12618498_R2.fastq.gz
SRR12618499_R1.fastq.gz
SRR12618499_R2.fastq.gz
SRR12618500_R1.fastq.gz
SRR12618500_R2.fastq.gz
SRR12618501_R1.fastq.gz
SRR12618501_R2.fastq.gz
```

```
sample,platform,fastq_1,fastq_2,lr_bam_file,bedfile
SRR12618483,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618483_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618483_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618484,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618484_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618484_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618485,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618485_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618485_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618486,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618486_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618486_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618487,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618487_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618487_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618488,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618488_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618488_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618489,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618489_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618489_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618490,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618490_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618490_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618491,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618491_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618491_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618492,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618492_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618492_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618493,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618493_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618493_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618494,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618494_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618494_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618495,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618495_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618495_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618496,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618496_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618496_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618497,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618497_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618497_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618498,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618498_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618498_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618499,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618499_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618499_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618500,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618500_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618500_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
SRR12618501,illumina,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618501_R1.fastq.gz,/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/SRR12618501_R2.fastq.gz,,https://raw.githubusercontent.com/artic-network/primer-schemes/master/nCoV-2019/V3/nCoV-2019.primer.bed
```

# Application Cont.

```
#!/bin/bash
#SBATCH --account=bphl-umbrella
#SBATCH --qos=bphl-umbrella
#SBATCH --job-name=aquascope
#SBATCH --cpus-per-task=8
#SBATCH --mem=64gb
#SBATCH --time=48:00:00
#SBATCH --output=aquascope.%j.out
#SBATCH --error=aquascope.%j.err
#SBATCH --mail-user=nikhil.yengala@flhealth.gov
#SBATCH --mail-type=FAIL,END

module purge
module load apptainer
module load conda

# Activate conda env with Nextflow
conda activate Nextflow

BASE="/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013"
samplesheet="${BASE}/input/samplesheet.csv"
outdir="${BASE}/results"
workdir="${outdir}/work"

if [[ ! -f "$samplesheet" ]]; then
    echo "ERROR: samplesheet not found: $samplesheet"
    exit 1
fi

mkdir -p "$outdir" "$workdir"

export APPTAINER_CACHEDIR="${BASE}/.apptainer_cache"
export NXF_APPTAINER_CACHEDIR="${BASE}/.apptainer_cache"

echo "Launching Aquascope..."
nextflow run main.nf \
    -entry AQUASCOPE \
    -profile singularity \
    --input "$samplesheet" \
    --outdir "$outdir" \
    -work-dir "$workdir" \
    -resume
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/
Name
├── ..
├── assets
├── bin
├── conf
├── docs
├── input
├── misc
├── modules
├── results
├── subworkflows
├── test
├── workflows
├── aquascope.14354509.err
├── aquascope.14354509.out
├── CHANGELOG.md
├── CITATIONS.md
├── code-of-conduct.md
├── CONTRIBUTING.md
├── DISCLAIMER.md
├── LICENSE
├── main.nf
├── mkdocs.yml
├── modules.json
├── nextflow.config
├── nextflow_schema.json
├── open_practices.md
├── README.md
├── ro-crate-metadata.json
├── rules_of_behavior.md
├── run_workflow.sh
├── tower.yml
└── VERSION
```

# Application Cont.

```
/.../n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/ALIGNMENT/
Name
└─ ..
   ├── SRR12618483.bam
   ├── SRR12618484.bam
   ├── SRR12618485.bam
   ├── SRR12618486.bam
   ├── SRR12618487.bam
   ├── SRR12618488.bam
   ├── SRR12618489.bam
   ├── SRR12618490.bam
   ├── SRR12618491.bam
   ├── SRR12618492.bam
   ├── SRR12618493.bam
   ├── SRR12618494.bam
   ├── SRR12618495.bam
   ├── SRR12618496.bam
   ├── SRR12618497.bam
   ├── SRR12618498.bam
   ├── SRR12618499.bam
   ├── SRR12618500.bam
   └── SRR12618501.bam
```









**sbatch** run\_workflow.sh

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/
Name
└─ ..
   ├── assets
   ├── bin
   ├── conf
   ├── docs
   ├── input
   ├── misc
   ├── modules
   └── results
       ├── subworkflows
       ├── test
       ├── workflows
       ├── aquascope.14354509.err
       ├── aquascope.14354509.out
       ├── CHANGELOG.md
       ├── CITATIONS.md
       ├── code-of-conduct.md
       ├── CONTRIBUTING.md
       ├── DISCLAIMER.md
       ├── LICENSE
       ├── main.nf
       ├── mkdocs.yml
       ├── modules.json
       ├── nextflow.config
       ├── nextflow_schema.json
       ├── open_practices.md
       ├── README.md
       ├── ro-crate-metadata.json
       ├── rules_of_behavior.md
       ├── run_workflow.sh
       ├── tower.yml
       └── VERSION
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/
Name
└─ ..
   ├── ALIGNMENT
   ├── BAMQC
   ├── FASTP
   ├── FASTQC
   ├── MULTIQC
   ├── pipeline_info
   ├── PRIMER_TRIMMING
   ├── SAMTOOLS
   ├── TRIMMED_SORTED_BAM
   ├── VARIANTS
   └── work
```



\_\_\_\_\_



\_\_\_\_\_

Summary

<b>General</b>	
fastp version:	0.23.4 ( <a href="https://github.com/OpenGenome/fastp">https://github.com/OpenGenome/fastp</a> )
sequencing:	paired end (101 cycles + 101 cycles)
mean length before filtering:	101bp, 101bp
mean length after filtering:	70bp, 70bp
duplication rate:	48.3422284
Insert size peak:	32
Detected read1 adapter:	AGATCGGAAGAGCACACGTCTGAACTCCAGTCA
Detected read2 adapter:	AGATCGGAAGAGCGTCGTGTAGAGAAAGAAGTT

**Before filtering**

total reads:	2.443002 M
total bases:	246.743202 M
Q30 bases:	235.073983 M (95.270703%)
Q30 bases:	226.720291 M (91.885122%)
GC content:	52.929515%

**After filtering**

**Filtering result**

reads passed filters:	2.034962 M (83.297599%)
reads with low quality:	1.249000 K (0.051055%)
reads with too many N:	274 (0.011216%)
reads too short:	404.518000 M (16.440101%)




**Adapters**

**Adapter or bad ligation of read1**





Sequence	Occurrences
AGATCGGAAGAGCACACGTCTGAACTCCAGTCACTTCGATTCGTTATGCGC	19700
AGATCGGAAGAGCGTCGAACTCCAGTCACTTCGATTCGTTATGCGCTCTCTCT	16245
AGATCGGAAGAGCACACGTCTGAACTCCAGTCACTTCGATTCGTTATGCGCTCTCTCTGA	7657
AGATCGGAAGAGCACACGTCTGAACTCCAGTCACTTCGATTCGTTATGCGCTCTCTCTGA	360931
AGATCGGAAGAGCACACGTCTGAACTCCAGTCACTTCGATTCGTTATGCGCTCTCTCTGTA	100945
AGATCGGAAGAGCGTCGAACTCCAGTCACTTCGATTCGTTATGCGCTCTCTCTGTA	11100
AGATCGGAAGAGCGTCGAACTCCAGTCACTTCGATTCGTTATGCGCTCTCTCTGTA	119919
other adapter sequences	128640

**Adapter or bad ligation of read2**

Sequence	Occurrences
AGATCGGAAGAGCGTCGTGTAGAGAAAGATTCTGAGAGCGTTAGATCTCGTGTGCTCCGAGATCAT	21806
AGATCGGAAGAGCGTCGTGTAGAGAAAGATTCTGAGAGCGTTAGATCTCGTGTGCTCCGATCA	10403
AGATCGGAAGAGCGTCGTGTAGAGAAAGATTCTGAGAGCGTTAGATCTCGTGTGCTCCGATCA	121425
AGATCGGAAGAGCGTCGTGTAGAGAAAGATTCTGAGAGCGTTAGATCTCGTGTGCTCCGATCA	9384
AGATCGGAAGAGCGTCGTGTAGAGAAAGATTCTGAGAGCGTTAGATCTCGTGTGCTCCGATCA	8872
AGATCGGAAGAGCGTCGTGTAGAGAAAGATTCTGAGAGCGTTAGATCTCGTGTGCTCCGATCA	14454
other adapter sequences	479150



# Application Cont.



# Application Cont.



A modular tool to aggregate results from bioinformatics analyses across many samples into a single report.

This report has been generated by the `nf-core/aquascope` analysis pipeline. For information about how to interpret these results, please see the [documentation](#).

Report generated on 2025-09-30, 16:00 EDT based on data in: `/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/work/44/500341606450621906w2147470806`

Welcome! Not sure where to start? [Watch a tutorial video](#) (2:00)

## General Statistics


Sample Name	% GC	Ins. size	± 30X	Median cov	Mean cov	% On target	% Aligned	Top lineage	Top lineage %	% Duplication	Reads After Filtering	GC content	% PF	% Adapter	Dups	GC	Median len	Segs
SRR12618483_1	40%	215	93.6%	149 X	182.8 X	100.0%	2.7%	Other	100.0%	88.3%	2.0M	54.9%	83.3%	54.5%	92.7%	54.0%	100bp	1.0M
SRR12618483_2	40%	215	89.2%	97 X	119.3 X	100.0%	2.9%	Other	100.0%	84.0%	1.3M	53.1%	79.4%	69.2%	91.9%	54.0%	100bp	1.0M
SRR12618484_1	40%	215	96.3%	92 X	117.4 X	100.0%	1.1%	Other	100.0%	76.6%	3.1M	56.2%	87.3%	38.3%	94.0%	53.0%	34bp	0.6M
SRR12618484_2	40%	215	90.0%	102 X	111.4 X	100.0%	1.9%	Other	100.0%	80.5%	1.8M	56.2%	82.4%	50.1%	94.1%	51.0%	101bp	0.9M
SRR12618485_1	40%	215	96.3%	92 X	117.4 X	100.0%	1.1%	Other	100.0%	76.6%	3.1M	56.2%	87.3%	38.3%	91.6%	55.0%	101bp	1.8M
SRR12618485_2	40%	217	90.0%	102 X	111.4 X	100.0%	1.9%	Other	100.0%	80.5%	1.8M	56.2%	82.4%	50.1%	90.5%	55.0%	101bp	1.8M
SRR12618486_1	40%	217	90.0%	102 X	111.4 X	100.0%	1.9%	Other	100.0%	80.5%	1.8M	56.2%	82.4%	50.1%	93.3%	52.0%	101bp	1.1M
SRR12618486_2	40%	37	8.0%	0 X	9.9 X	100.0%	0.3%	A	84.6%	81.8%	1.7M	53.4%	82.3%	85.6%	91.4%	56.0%	100bp	0.9M
SRR12618487_1	40%	37	8.0%	0 X	9.9 X	100.0%	0.3%	A	84.6%	81.8%	1.7M	53.4%	82.3%	85.6%	84.2%	60.0%	101bp	1.0M
SRR12618487_2	40%	37	8.0%	0 X	9.9 X	100.0%	0.3%	A	84.6%	81.8%	1.7M	53.4%	82.3%	85.6%	92.4%	53.0%	32bp	0.9M

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/  
Name  
..  
ALIGNMENT  
BAMQC  
FASTP  
FASTQC  
MULTIQC  
pipeline_info  
PRIMER_TRIMMING  
SAMTOOLS  
TRIMMED_SORTED_BAM  
VARIANTS  
work
```

```
../n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/MULTIQC/  
Name  
..  
multiqc_data  
multiqc_plots  
multiqc_report.html  
versions.yml
```


## Cumulative genome coverage

Percentage of the reference genome with at least the given depth of coverage.



## GC content distribution


Each solid line represents the distribution of GC content of mapped reads for a given sample.



## Filtered Reads

Filtering statistics of sampled reads.

Percentages



# Application Cont.

## Nextflow workflow report

[sharp\_wozniak] *(resumed run)*

Workflow execution completed successfully!

**Run times**  
30-Sep-2025 15:40:56 - 30-Sep-2025 16:00:23 (duration: **19m 27s**)

251 succeeded

**Nextflow command**  
nextflow run main.nf -entry AQUASCOPE -profile singularity --input /blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/input/samplesheet.csv --outdir /blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results -work-dir /blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/work -resume

CPU-Hours	1.6
Launch directory	/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013
Work directory	/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/work
Project directory	/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013
Script name	main.nf
Script ID	267875cc96f133a8a8cf602aad079305
Workflow session	2f68cb22-1045-405e-bdeb-4df50a53c2b3
Workflow profile	singularity
Nextflow version	version 24.10.2, build 5932 (27-11-2024 21:23 UTC)

/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/

Name

..

ALIGNMENT

BAMQC

FASTP

FASTQC

MULTIQC

pipeline\_info

PRIMER\_TRIMMING

SAMTOOLS

TRIMMED\_SORTED\_BAM

VARIANTS

work

/.../n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/pipeline\_info/

Name

..

execution\_report\_2025-09-30\_15-40-55.html

execution\_timeline\_2025-09-30\_15-40-55.html

execution\_trace\_2025-09-30\_15-40-55.txt

params\_2025-09-30\_15-40-59.json

pipeline\_dag\_2025-09-30\_15-40-55.html

samplesheet.valid.csv

versions.yml




# Application Cont.

```
./repos/bphl-molecular/analysis/aquascope1013/results/PRIMER_TRIMMING/IVAR/
Name
..
SRR12618483.IVAR_trim.bam
SRR12618483.IVAR_trim.ivar.log
SRR12618484.IVAR_trim.bam
SRR12618484.IVAR_trim.ivar.log
SRR12618485.IVAR_trim.bam
SRR12618485.IVAR_trim.ivar.log
SRR12618486.IVAR_trim.bam
SRR12618486.IVAR_trim.ivar.log
SRR12618487.IVAR_trim.bam
SRR12618487.IVAR_trim.ivar.log
SRR12618488.IVAR_trim.bam
SRR12618488.IVAR_trim.ivar.log
SRR12618489.IVAR_trim.bam
SRR12618489.IVAR_trim.ivar.log
SRR12618490.IVAR_trim.bam
SRR12618490.IVAR_trim.ivar.log
SRR12618491.IVAR_trim.bam
SRR12618491.IVAR_trim.ivar.log
SRR12618492.IVAR_trim.bam
SRR12618492.IVAR_trim.ivar.log
SRR12618493.IVAR_trim.bam
SRR12618493.IVAR_trim.ivar.log
SRR12618494.IVAR_trim.bam
SRR12618494.IVAR_trim.ivar.log
SRR12618495.IVAR_trim.bam
SRR12618495.IVAR_trim.ivar.log
SRR12618496.IVAR_trim.bam
SRR12618496.IVAR_trim.ivar.log
SRR12618497.IVAR_trim.bam
SRR12618497.IVAR_trim.ivar.log
SRR12618498.IVAR_trim.bam
SRR12618498.IVAR_trim.ivar.log
SRR12618499.IVAR_trim.bam
SRR12618499.IVAR_trim.ivar.log
SRR12618500.IVAR_trim.bam
SRR12618500.IVAR_trim.ivar.log
SRR12618501.IVAR_trim.bam
SRR12618501.IVAR_trim.ivar.log
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/
Name
..
ALIGNMENT
BAMQC
FASTP
FASTQC
MULTIQC
pipeline_info
PRIMER_TRIMMING
SAMTOOLS
TRIMMED_SORTED_BAM
VARIANTS
work
```

```
./repos/bphl-molecular/analysis/aquascope1013/results/PRIMER_TRIMMING/
Name
..
IVAR
```



# Application Cont.

```
54907 + 0 in total (QC-passed reads + QC-failed reads)
54880 + 0 primary
0 + 0 secondary
27 + 0 supplementary
0 + 0 duplicates
0 + 0 primary duplicates
54907 + 0 mapped (100.00% : N/A)
54880 + 0 primary mapped (100.00% : N/A)
54880 + 0 paired in sequencing
27511 + 0 read1
27369 + 0 read2
54697 + 0 properly paired (99.67% : N/A)
54745 + 0 with itself and mate mapped
135 + 0 singletons (0.25% : N/A)
0 + 0 with mate mapped to a different chr
0 + 0 with mate mapped to a different chr (mapQ>=5)
```

```
# This file was produced by samtools stats (1.21+htslib-1.21) and can be plotted using plot-bamstats
# This file contains statistics for all reads.
# The command line was: stats --threads 1 --reference SARS-CoV-2.reference.fasta SRR12618483.trimmed.sorted.bam
# CHK, Checksum [2]Read Names [3]Sequences [4]Qualities
# CHK, CRC32 of reads which passed filtering followed by addition (32bit overflow)
CHK 1102a2d3 171eb6fb 228185ca
# Summary Numbers. Use 'grep ^SN | cut -f 2-' to extract this part.
SN raw total sequences: 54880 # excluding supplementary and secondary reads
SN filtered sequences: 0
SN sequences: 54880
SN is sorted: 1
SN 1st fragments: 27511
SN last fragments: 27369
SN reads mapped: 54880
SN reads mapped and paired: 54745 # paired-end technology bit set + both mates mapped
SN reads unmapped: 0
SN reads properly paired: 54697 # proper-pair bit set
SN reads paired: 54880 # paired-end technology bit set
SN reads duplicated: 0 # PCR or optical duplicate bit set
SN reads MQ0: 1 # mapped and MQ=0
SN reads QC failed: 0
SN non-primary alignments: 0
SN supplementary alignments: 27
SN total length: 5462200 # ignores clipping
SN total first fragment length: 2738133 # ignores clipping
SN total last fragment length: 2724067 # ignores clipping
SN bases mapped: 5462200 # ignores clipping
SN bases mapped (cigar): 5312914 # more accurate
SN bases trimmed: 0
SN bases duplicated: 0
SN mismatches: 9836 # from NM fields
SN error rate: 1.851338e-03 # mismatches / bases mapped (cigar)
SN average length: 100
SN average first fragment length: 100
SN average last fragment length: 100
SN maximum length: 101
SN maximum first fragment length: 101
SN maximum last fragment length: 101
SN average quality: 39.0
SN insert size average: 202.0
SN insert size standard deviation: 69.2
SN inward oriented pairs: 26891
SN outward oriented pairs: 432
SN pairs with other orientation: 5
SN pairs on different chromosomes: 0
SN percentage of properly paired reads (%): 99.7
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/
Name
├── ..
├── ALIGNMENT
├── BAMQC
├── FASTP
├── FASTQC
├── MULTIQC
├── pipeline_info
├── PRIMER_TRIMMING
├── SAMTOOLS
├── TRIMMED_SORTED_BAM
├── VARIANTS
└── work
```

```
/.../n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/SAMTOOLS/
Name
├── ..
├── SRR12618483.flagstat
├── SRR12618483.idxstats
├── SRR12618483.stats
├── SRR12618484.flagstat
├── SRR12618484.idxstats
└── SRR12618484.stats
```

# Application Cont.

```
./../repos/bphl-molecular/analysis/aquascope1013/results/TRIMMED_SORTED_BAM/  
Name  
..  
SRR12618483.trimmed.sorted.bam  
SRR12618483.trimmed.sorted.bam.bai  
SRR12618484.trimmed.sorted.bam  
SRR12618484.trimmed.sorted.bam.bai  
SRR12618485.trimmed.sorted.bam  
SRR12618485.trimmed.sorted.bam.bai  
SRR12618486.trimmed.sorted.bam  
SRR12618486.trimmed.sorted.bam.bai  
SRR12618487.trimmed.sorted.bam  
SRR12618487.trimmed.sorted.bam.bai  
SRR12618488.trimmed.sorted.bam  
SRR12618488.trimmed.sorted.bam.bai  
SRR12618489.trimmed.sorted.bam  
SRR12618489.trimmed.sorted.bam.bai  
SRR12618490.trimmed.sorted.bam  
SRR12618490.trimmed.sorted.bam.bai  
SRR12618491.trimmed.sorted.bam  
SRR12618491.trimmed.sorted.bam.bai  
SRR12618492.trimmed.sorted.bam  
SRR12618492.trimmed.sorted.bam.bai  
SRR12618493.trimmed.sorted.bam  
SRR12618493.trimmed.sorted.bam.bai  
SRR12618494.trimmed.sorted.bam  
SRR12618494.trimmed.sorted.bam.bai  
SRR12618495.trimmed.sorted.bam  
SRR12618495.trimmed.sorted.bam.bai  
SRR12618496.trimmed.sorted.bam  
SRR12618496.trimmed.sorted.bam.bai  
SRR12618497.trimmed.sorted.bam  
SRR12618497.trimmed.sorted.bam.bai  
SRR12618498.trimmed.sorted.bam  
SRR12618498.trimmed.sorted.bam.bai  
SRR12618499.trimmed.sorted.bam  
SRR12618499.trimmed.sorted.bam.bai  
SRR12618500.trimmed.sorted.bam  
SRR12618500.trimmed.sorted.bam.bai  
SRR12618501.trimmed.sorted.bam  
SRR12618501.trimmed.sorted.bam.bai
```

```
/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/  
Name  
..  
ALIGNMENT  
BAMQC  
FASTP  
FASTQC  
MULTIQC  
pipeline_info  
PRIMER_TRIMMING  
SAMTOOLS  
TRIMMED_SORTED_BAM  
VARIANTS  
work
```





# Application Cont.

	SRR12618483.variants.tsv
summarized	[['Other', 1.0000000012495296]]
lineages	B.6 B.45 B.1.301 B.1.369 B.1.1.228 B.1.578
abundances	0.28440400 0.20719400 0.20563400 0.13082300 0.12865500 0.04329000
resid	4.423475409
coverage	97.40828679

REGION	POS	REF	ALT	REF_DP	REF_RV	REF_QUAL	ALT_DP	ALT_RV	ALT_QUAL	ALT_FREQ	TOTAL_DP	PVAL	PASS	GFF_FEATURE	REF_CODON	REF_AA	ALT_CODON	ALT_AA	POS_AA
MN908947.3	213	G	A	53	0	39	1	0	40	0.0185185	54	0.504673	FALSE	NA	NA	NA	NA	NA	NA
MN908947.3	241	C	T	80	0	39	54	24	69	0.402985	134	8.20E-20	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	329	C	T	79	6	43	1	1	78	0.0125	80	0.503145	FALSE	NA	NA	NA	NA	NA	NA
MN908947.3	437	A	C	165	123	39	9	0	39	0.0517241	174	0.00180324	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	442	C	A	135	85	39	38	38	39	0.219653	173	4.16E-13	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	662	G	A	181	56	39	2	1	34	0.010929	183	0.250685	FALSE	NA	NA	NA	NA	NA	NA
MN908947.3	746	T	#NAME?	349	166	47	19	0	20	0.0544413	349	0.000568383	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	748	G	A	329	147	48	21	20	37	0.06	350	3.61E-07	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	752	A	-C	346	166	47	52	0	20	0.149856	347	4.90E-12	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	906	C	T	532	420	47	17	9	39	0.0309654	549	6.83E-06	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	942	G	A	296	180	52	24	22	45	0.075	320	3.96E-08	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	1018	G	C	114	8	42	9	0	39	0.0731707	123	0.0017421	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	1059	C	T	33	0	39	14	14	38	0.297872	47	2.30E-05	TRUE	NA	NA	NA	NA	NA	NA
MN908947.3	1212	G	A	254	137	44	31	31	39	0.108772	285	2.08E-10	TRUE	NA	NA	NA	NA	NA	NA

MN908947.3	1	A	0
MN908947.3	2	T	0
MN908947.3	3	T	0
MN908947.3	4	A	0
MN908947.3	5	A	0
MN908947.3	6	A	0
MN908947.3	7	G	0
MN908947.3	8	G	0
MN908947.3	9	T	0
MN908947.3	10	T	0
MN908947.3	11	T	0
MN908947.3	12	A	0
MN908947.3	13	T	0
MN908947.3	14	A	0

/blue/bphl-florida/n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/	
Name	^
..	
ALIGNMENT	
BAMQC	
FASTP	
FASTQC	
MULTIQC	
pipeline_info	
PRIMER_TRIMMING	
SAMTOOLS	
TRIMMED_SORTED_BAM	
VARIANTS	
work	

/.../n.yengalareddy/repos/bphl-molecular/analysis/aquascope1013/results/VARIANTS/	
Name	^
..	
FREYJA	

/.../repos/bphl-molecular/analysis/aquascope1013/results/VARIANTS/FREYJA/	
Name	^
..	
Demix	
freyja_db	
VarCalls	
versions.yml	

/.../repos/bphl-molecular/analysis/aquascope1013/results/VARIANTS/FREYJA/Demix/	
Name	^
..	
SRR12618483.tsv	
SRR12618484.tsv	
SRR12618485.tsv	

/.../repos/bphl-molecular/analysis/aquascope1013/results/VARIANTS/FREYJA/freyja_db/	
Name	^
..	
curated_lineages.json	
lineages.yml	
usher_barcode.feather	

/.../repos/bphl-molecular/analysis/aquascope1013/results/VARIANTS/FREYJA/VarCalls/	
Name	^
..	
SRR12618483.depth.tsv	
SRR12618483.variants.tsv	
SRR12618484.depth.tsv	

# Conclusion



Fundamentals of  
Aquascope



Installation and setup  
of Aquascope in HPG



Successfully executed  
job query for  
Aquascope



Generated output files



# Advanced Molecular Detection

## Southeast Region Bioinformatics

# Questions?

[bphl-sebioinformatics@flhealth.gov](mailto:bphl-sebioinformatics@flhealth.gov)

**Molly Mitchell, PhD**

Bioinformatics Supervisor

[Molly.Mitchell@flhealth.gov](mailto:Molly.Mitchell@flhealth.gov)

**Nikhil Reddy, MS**

Bioinformatician

[Nikhil.Yengala@flhealth.gov](mailto:Nikhil.Yengala@flhealth.gov)

**Sam Bernhoft, MPH**

Bioinformatician

[Samantha.bernhoft@flhealth.gov](mailto:Samantha.bernhoft@flhealth.gov)