SeqSender

Advanced Molecular Detection

Southeast Region Bioinformatics

This resource was made possible through funding provided under the Epidemiology and Laboratory Capacity for Prevention and Control of Emerging Infectious Diseases (ELC) Cooperative Agreement (CK24-0002), Project D: Advanced Molecular Detection to the Florida Department of Health. The conclusions, findings, and opinions expressed by authors do not necessarily reflect the official position of the U.S. Department of Health and Human Services, the Public Health Service, or the Centers for Disease Control and Prevention

OVERVIEW

<u>Purpose</u>

• A CLI tool developed by the CDC to streamline and automate the submission of viral and pathogen WGS data to public databases such as GenBank, SRA, GISAID and BioSample.

<u>Usage</u>


• Used by public health labs, researchers, and sequencing centers to generate standardized submission files and upload assembled genomes with metadata to repositories.

Dependencies

Conda



WORKFLOW






Application

Objective

Use a SARS-CoV-2 dataset and generate an output from the SeqSender pipeline






```
cd /blue/bphl-<state>/<user>/repos/bphl-
molecular/
git clone https://github.com/CDCgov/seqsender
mkdir analysis/
cd analysis/
cp /blue/bphl-<state>/<user>/repos/bphl-
molecular/ seqsender/
```



```
GNU nano 5.6.1
                                                       run_segsender.sh
#!/usr/bin/bash
#SBATCH --account=bphl-umbrella
#SBATCH --gos=bphl-umbrella
#SBATCH --job-name=segsender
#SBATCH --ntasks=1
#SBATCH --cpus-per-task=20
#SBATCH --mem=20gb
#SBATCH --time=48:00:00
#SBATCH --output=segsender.%j.out
#SBATCH --error=segsender.%j.err
python seqsender.py submit
 --organism COV \
  --submission name run1
  --submission_dir output
  --config_file test_data/COV/cov_config.yaml \
  --metadata_file test_data/COV/cov_genbank_metadata.csv \
  --fasta_file test_data/COV/cov_sequence.fasta \
```


Submitting To NCB

Enter your username and password in the config file







conda env –f env.yaml
conda activate seqsender
shatch run_seqsender.sh






















Conclusion







Advanced Molecular Detection Southeast Region Bioinformatics

Questions?

bphl-sebioinformatics@flhealth.gov

Molly Mitchell, PhD
Bioinformatics Supervisor
Molly.Mitchell@flhealth.gov

Nikhil Reddy, MS
Bioinformatician
Nikhil.Yengala@flhealth.gov

Sam Bernhoft, MPH
Bioinformatician
Samantha.bernhoft@flhealth.gov